# RADIESSE® + LIDOCAINE

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---

### DESCRIPTION

RADIESSE® (+) Lidocaine injectable implant is a subdermal, semi-solid, cohesive and resorbable implant, non-pyrogenic, latex free, sterile. The principal component is synthetic calcium hydroxyapatite, a biocompatible material.

The semi-solid nature of the implant originates from the suspension of calcium hydroxyapatite in a gel-based vehicle, consisting mainly of water (sterile water for injection USP), sodium carboxymethylcellulose (USP), glycerin (USP) and lidocaine hydrochloride. The CaHA-containing gel is formulated to stimulate an increase in the production of collagen and elastin in the treated region. The elastic and collagen fibers are in the middle of the extracellular matrix (ECM) of the dermis and are responsible for important biomechanical properties of the skin. The result is a long-lasting restoration and volumization, with improvements in laxity and skin tightening in the injected region by collagen biostimulation, though not permanent.

RADIESSE® (+) Lidocaine injectable implant has particle size range of 25–45 microns and should be injected with a 27 gauge Inner Diameter (I.D.) needle. **The use of needles of caliber less than 27G may increase the incidence of occlusion.**

### PRESENTATIONS

RADIESSE® (+) Lidocaine is supplied in a sterilized filled syringe for single use. Each carton contains one syringe of 1.5cc, one instruction for use, 3 traceability labels, two sterile needles, for single use only. The dimensions of the needles are specified on the outer carton.

### COMPOSITION

The RADIESSE® (+) Lidocaine Injectable Implant consists of approximately 30% calcium hydroxyapatite per volume and 70% gel as vehicle, by volume.

Composition: calcium hydroxyapatite, sterile water for injection, glycerin, sodium carboxymethylcellulose, dibasic sodium phosphate, sodium phosphate monobasic as buffers and lidocaine hydrochloride (0.3%).

### INTENDED USE

The RADIESSE® (+) Lidocaine Injectable Implant is indicated for:

- Plastic and reconstructive surgery, including deep filling of the subdermal soft tissues, such as correction of wrinkles and moderate to severe furrows (such as nasolabial folds), contouring and improvements in laxity and skin tightening in the injected region by biostimulation of collagen.

It acts by adding the volume of tissue, restoring or increasing the contours of the skin, depending on the desired amount of correction.

The CaHA gel matrix is formulated to stimulate an increase in the production of collagen with improvements in laxity and skin tightening in the injected region (areas such as facial, abdomen, thighs and brachial zone, upper arms, neck, décolletage, etc.). Collagen and elastin are components of the extracellular matrix (ECM) that confer important biomechanical properties to the skin. The result is a long-lasting restoration and volumization, with improvements in laxity and skin tightening in the injected region by collagen biostimulation, though not permanent.

### CONTRAINDICATIONS

RADIESSE® (+) Lidocaine is contraindicated for patients with:

- Severe allergies manifested by a history of anaphylaxis, or history or presence of multiple severe allergies.
- Known hypersensitivity to lidocaine or anesthetics of the amide type.
- If the patient has active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- Known hypersensitivity to any of the components.
- In patients prone to developing inflammatory skin conditions or those patients with a tendency for developing hypertrophic scars and keloids.
- This product shall not be implanted in the epidermis or used as a skin replacement. Implantation into the epidermis or superficial dermis could lead to complications such as fistula formation, infections, extrusions, nodule formation and induration.
- This product is not intended to be used glabellar folds. A higher incidence of localized necrosis has been associated with glabellar injection. Complications associated with another injectables indicate that forceful injection into superficial dermal vessels of the glabellar area can result in vascular occlusion.

- In the presence of foreign bodies such as liquid silicone or other particulate materials.
- In areas where there is inadequate coverage of healthy, well vascularized tissue.
- In patients with systemic disorders which cause poor wound healing or will lead to tissue deterioration over the implant.
- Patients with bleeding disorders.

## **WARNINGS**

- Special attention must be taken in order to avoid the implant into the vasculature because it may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention.
- Use of RADIESSE® (+) Lidocaine injectable implant in any person with active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- Do not overcorrect (overfill) a contour deficiency because the depression should gradually improve within several weeks as the treatment effect of RADIESSE® (+) Lidocaine injectable implant occurs. Refer to the individualization of treatment sections for details.
- The safety and effectiveness for use in the lips has not been established. There have been published reports of nodules associated with the use of RADIESSE® (+) Lidocaine injectable implant injected into the lips.
- Injection procedure reactions have been observed consisting mainly of short-term (i.e., <7 days) such as bruising, redness and swelling. Refer to the Adverse Events sections for details.

## **PRECAUTIONS**

- In order to minimize the risks of potential complications, the product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. Health care practitioners should fully familiarize themselves with the product, the product educational materials and the entire package insert. Health care practitioners are encouraged to discuss all potential risks of soft tissue injections with their patients prior to treatment and ensure that the patients are aware of signs and symptoms of potential complications.
- Universal precautions must be observed when there is a potential for contact with patient body fluids. The injection session must be conducted with aseptic technique.
- RADIESSE® (+) Lidocaine injectable implant is packaged for single patient use. Do not resterilize. Do not use if package is opened or damaged. Do not use if the syringe end cap or syringe plunger is not in place.
- The calcium hydroxyapatite (CaHA) particles of RADIESSE® (+) Lidocaine injectable implant are radiopaque and are clearly visible on CT Scans and may be visible in standard, plain radiography. Patients need to be informed of the radio-opaque nature of RADIESSE® (+) Lidocaine injectable implant, so that they can inform their primary care health professionals as well as radiologists. In a radiographic study of 58 patients, there was no indication of RADIESSE® injectable implant potentially masking abnormal tissues or being interpreted as tumors in CT Scans.
- The long-term safety, more than 3 years, of RADIESSE® (+) Lidocaine injectable implant has not been evaluated in clinical trials.
- The safety of RADIESSE® (+) Lidocaine injectable implant in patients with increased susceptibility to keloid formation and hypertrophic scarring has not been studied.
- As with all transcutaneous procedures, RADIESSE® (+) Lidocaine injectable implant injection carries a risk of infection. Infection may necessitate attempted surgical removal of RADIESSE® (+) Lidocaine. Standard precautions associated with injectable materials should be followed.
- Safety of RADIESSE® (+) Lidocaine injectable implant for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established.
- The patient should be informed that he or she should minimize exposure of the treated area to extensive sun or heat exposure for approximately 24 hours after treatment or until any initial swelling and redness has resolved.
- Patients who are using medications that can prolong bleeding, such as aspirin or warfarin, may, as with any injection, experience increased bruising or bleeding at the injection site.
- Safety and effectiveness in the periorbital area has not been established.
- No studies of interactions of RADIESSE® (+) Lidocaine injectable implant with drugs or other substances or implants have been conducted.

- In healthy adults, the maximum total dose of lidocaine (without epinephrine) should not exceed 300 mg (or 4.5 mg / kg) per session. Overdose with lidocaine usually results in changes in the central nervous system or cardiovascular toxicity. When used concomitantly (topical administration), the total dose of lidocaine administered should be taken into account. Concomitant use of other local anesthetic agents or structurally related to amide-type local anesthetics should also be considered given that toxic systemic effects may be additive.
- The safety of RADIESSE® (+) Lidocaine injectable implant with concomitant dermal therapies such as epilation, UV irradiation, or laser, mechanical or chemical peeling procedures has not been evaluated in controlled clinical trials.
- If laser treatment, chemical peeling, or any other procedure based on active dermal response is considered after treatment with RADIESSE® (+) Lidocaine injectable implant, there is a possible risk of eliciting an inflammatory reaction at the implant site. This also applies if RADIESSE® (+) injectable implant is administered before the skin has healed completely after such a procedure.
- Injection of RADIESSE® (+) Lidocaine injectable implant into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes.
- To help avoid needle breakage, do not attempt to straighten a bent needle. Discard it and complete the procedure with a replacement needle.
- After use, treatment syringes and needles may be potential biohazards. Handle accordingly and dispose of in accordance with accepted medical practice and applicable local, state and federal requirements.
- Care should be taken to determine the risk versus the benefit for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies, and with patients who are receiving concomitant treatment with methemoglobin-inducing agent.

## ADVERSE EVENTS

Adverse events observed in clinical trials with the RADIESSE® (+) Lidocaine Injection Implant include reactions due to injection, and to the needle.

In a randomized, controlled trial of the RADIESSE® (+) Lidocaine Injection Implant for the treatment of nasolabial folds, patients received RADIESSE® (+) Lidocaine on one side of the face / groove (treatment group) and on the other side RADIESSE® (control). Adverse events reported during this study were generally expected to be mild in nature and short in duration. The majority of adverse events were reported through patient diaries.

Table 1 summarizes the number of adverse events reported in the subject diaries. Swelling and redness were the most frequently reported adverse events. There was no significant difference in adverse event rates between the nasolabial folds with RADIESSE® (+) Lidocaine and the nasolabial folds with RADIESSE® (without lidocaine).

No vascular compromise events occurred in the RADIESSE® (+) Lidocaine injectable implant injections.

In the RADIESSE® (without lidocaine) injections, two (2/101, 2%) vascular compromise events occurred, requiring treatment to resolve.

**Table 1. Adverse Events Reported in Subject Diaries Over the 4-Weeks Study Period**  
N = 202 Folds

Event Type *	N (%)	
	RADIESSE® (+) N = 101 Folds	RADIESSE® N = 101 Folds
Bruising	44 (43.6%)	48 (47.5%)
Itching	37 (36.6%)	34 (33.7%)
Pain	48 (47.5%)	57 (56.4%)
Redness	66 (65.3%)	71 (70.3%)
Swelling	90 (89.1%)	92 (91.1%)
Blanching	5 (5.0%)	8 (7.9%)

\* "Other" adverse events reported by 19 subjects for both RADIESSE® injectable implant and Control include numbness, tenderness, lumps, bumps and discomfort.

Of the 13 blanching events described in Table 1, two were associated with previously described vascular compromise events. The remaining 11 were not determined to be vascular compromise events.

Table 2 summarizes the number of adverse events reported by the investigators. As with the patients, swelling and redness were the most frequently reported adverse events.

**Table 2. Adverse Events Reported by Investigators**  
N = 101 Subjects

Event Type *	N (%)	
	RADIESSE® (+)	RADIESSE®
Bruising	20 (19.8%)	18 (17.8%)
Swelling	58 (57.4%)	55 (54.5%)
Erosion	0	1 (1.0%)**
Redness	51 (50.5%)	50 (49.5%)
Infection	1 (1.0%)†	0
Needle Jamming	3 (3.0%)‡	0
Nodule	0	0
Pain	0	1 (1.0%)
Vascular Compromise	0	1 (1.0%)

\*"Other" adverse events for both RADIESSE® injectable implant and Control include tenderness and tingling sensation.

\*\* Associated with vascular compromise event.

† Herpes simplex infection unrelated to study device.

‡ Returned product evaluation regarding one subject concluded that needle jam might have resulted from needle incompletely screwed onto the RADIESSE® syringe's luer threads.

The use of lidocaine is well established as a local and regional anesthetic. Some potential side effects of lidocaine are dizziness, nervousness, apprehension, euphoria, confusion, loss of consciousness, drowsiness, ringing in the ears, blurred vision or double vision, vomiting, heat sensations, cold or numbness, spasms, tremors, seizures, respiratory depression and stopping, slow heartbeat, hypotension, and cardiovascular collapse (which can lead to cardiac arrest).

## POST MARKETING SURVEILLANCE

The following adverse events were received from post surveillance for the injectable RADIESSE® implant. They were not observed in clinical trials. They are voluntary reports so it is not possible to estimate their frequency or to establish a causal relationship with RADIESSE® (+) Lidocaine: pruritus, itching, angioedema, inflammation, necrosis, granuloma, nodules, induration, erythema, skin discoloration, pustule, skin pallor, hair loss, paresthesia, ptosis, pain, headache, swelling, asymmetry, abscess, herpetic infection including herpes simplex and herpes zoster, hematoma, blanching, blistering, dizziness, festoons, flu-like symptoms, Guillain-Barre syndrome, tachypnea, ischemic reaction, lymphoid hyperplasia, nausea, pericarditis, scarring, sensitivity to cold, vascular occlusion/obstruction, vascular compromise, ocular ischemia, diplopia, visual impairment/blindness, facial muscle paralysis, Bell's palsy.

The following interventions have been reported: antibiotics, anti-inflammatories, corticosteroids, anti-histamines, analgesics, massage, warm compress, excision, drainage, and surgery. This information does not constitute and is not intended to be medical advice, a recommendation on how to treat an adverse event or an exhaustive list of possible interventions. Health professionals should evaluate each case on an individual basis, and independently determine, based on their professional experience, what treatment(s) are appropriate, if any, for their patients.

## INDIVIDUALIZATION OF TREATMENT

Before treatment, the patient's suitability for the treatment and the patient's need for pain relief should be assessed. The outcome of treatment with RADIESSE® (+) Lidocaine injectable implant will vary between patients. In some instances, additional treatments may be necessary depending on the size of the defect and the needs of the patient.

## DIRECTIONS FOR USE

### GENERAL

The following is required for the percutaneous injection procedure:

- RADIESSE® (+) Lidocaine injectable implant syringe(s)
  - 25 gauge OD - 27 gauge ID needle(s) with Luer lock fittings
1. Prepare patient for percutaneous injection using standard methods. The treatment injection site should be marked and prepared with a suitable antiseptic for each patient.
  2. Prepare the syringes of RADIESSE® (+) Lidocaine injectable implant and the injection needle(s) before the percutaneous injection. A new injection needle may be used for each syringe, or the same injection needle may be connected to each new syringe for each patient.
  3. Remove foil pouch from the carton. Open the foil pouch by tearing at the notches (marked 1 and 2), and remove the syringe. There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.

4. Peel or twist apart the needle packaging to expose the hub.
5. Remove the Luer syringe cap from the distal end of the syringe prior to attaching the needle. The syringe of RADIESSE® (+) Lidocaine injectable implant can then be twisted onto the Luer lock fitting of the needle taking care not to contaminate the needle. Discard needle package. The needle must be tightened securely to the syringe. If excess implant is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until RADIESSE® (+) Lidocaine injectable implant extrudes from the end of the needle. If leakage is noted at the Luer fitting, it may be necessary to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the needle.
6. Locate the initial site for the implant. Scar tissue and cartilage may be difficult or impossible to treat. Avoid passing through these tissue types when advancing the injection needle during the procedure.
7. The amount injected will vary depending on the site and extent of the restoration or filling desired. RADIESSE® (+) Lidocaine injectable implant should be injected subdermally.
8. Use a 1:1 correction factor. No overcorrection is needed.
9. Insert needle with bevel down at approximately a 30° angle to the skin. Needle should slide under the dermis to the point you wish to begin the injection. This should be easily palpable with the non-dominant hand.
10. If significant resistance is encountered when pushing the plunger, the injection needle may be moved slightly to allow easier placement of the material or it may be necessary to change the injection needle. Needle jams are more likely with use of needles smaller than 27 gauge ID.
11. Advance the needle into the subdermis to the starting location. Carefully push the plunger of the RADIESSE® (+) Lidocaine injectable implant syringe to start the injection and slowly inject the implant material in linear threads while withdrawing the needle. Continue placing additional lines of material until the desired level of correction is achieved.
12. Apply slow continuous even pressure to the syringe plunger to inject the implant as you withdraw the needle. The implant material should be completely surrounded by soft tissue without leaving globular deposits. The injected area may be massaged as needed to achieve even distribution of the implant.
13. Use once and discard in accordance with local safety standards.

**Note: Health professionals should evaluate each case, individually, and independently determine, based on their professional experience, the appropriate treatment (s), if any, for their patients. Thus, the dilution of CaHA to the treatment may change according to the area treated, the degree of laxity and the thickness of the skin.**

#### **INSTRUCTIONS FOR TRACEABILITY LABELS**

A set of labels is supplied inside the package. The labels on the aluminum packaging are used to ensure product identification and traceability. These labels should be affixed as follows: one in the medical record, another in the document to be delivered to the patient, and the third in the control of the clinic / hospital.

#### **PATIENT COUNSELING INFORMATION**

For about 24 hours after injection, the following should be avoided in the treated area;

- Significant movements or massages in the treated area
- Applying makeup
- Extensive exposure to sun or heat

In addition, the Health professional's instructions regarding care of the injection site should be followed, and may include the use of cold packs and / or medications to reduce swelling and bruising.

#### **GENERAL INFORMATION**

##### **SUPPLIER**

The RADIESSE® (+) Lidocaine Injectable Implant is supplied sterile and non-pyrogenic in a syringe packaged in an aluminum foil pack, placed inside a carton for easy storage. Each unit includes a 1.5cc filled syringe of RADIESSE® (+) Lidocaine Injectable Implant. The degree of accuracy of syringe graduations is  $\pm 0.025\text{cc}$ . Do not use the product if the package and / or syringe is damaged or if the protective cap or plunger of the syringe is not intact.

**The contents of the syringe are intended for only one patient in a single treatment and cannot be re-sterilized. Re-use can compromise the functional properties of the device and / or lead to device failure. Re-use can also create a risk of contamination of the device and / or cause infection or cross-infection including, but not limited to, transmission of infectious diseases and transfer of blood between patients. All this can, in turn, lead to injury, illness or death of the patient.**

##### **STORAGE**

RADIESSE® (+) Lidocaine injectable implant should be stored at a controlled room temperature between 15°C and 32°C (59°F and 90°F). Do not use if the expiration date has been exceeded. The expiration date is printed on the product label.

## DISPOSAL

Used and partially used syringes and injection needles could be biohazardous and should be handled and disposed of in accordance with facility medical practices and local, state or federal regulations.

## WARRANTY

Merz North America, Inc. warrants that reasonable care has been exercised in the design and manufacture of this product.

**THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ITS PARTICULAR PURPOSE.**

Handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Merz North America, Inc.'s control directly affect the product and the results obtained from its use. Merz North America, Inc.'s obligation under this warranty is limited to the replacement of this product and Merz North America, Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly, arising from the use of this product. Merz North America, Inc. neither assumes, nor authorizes any person to assume for Merz North America, Inc., any other or additional liability or responsibility in connection with this product.

## SYMBOLS

 Caution.	 Consult instructions for use.
 Do not use if package is damaged.	 Single use product. Do not re-use. Forbidden reprocess.
 Not Made with Natural Rubber Latex.	 Use by date.
 Date of manufacture.	 Manufacturer.
 Lot.	 Sterile. Sterile using steam.
 15°C (59°F) – 32°C (90°F) Temperature limit of storage: 15°C – 32°C	

## Manufacturer

Merz North America, Inc.  
4133 Courtney St., Suite 10 - Franksville, WI 53126 USA

Country of Origin: United States of America.

Registered/ Imported and Distributed by: Merz Farmacêutica Comercial Ltda.

Headquarters: CNPJ: 11.681.446/0001-45 Avenida Doutor Cardoso de Melo nº 1184 Conj.131 e 132, Vila Olímpia, São Paulo/SP CEP 04548-004

Branches: CNPJ: 11.681.446/0003-07

CNPJ: 11.681.446/0005-79

Technical Responsible: Angelica Santinelli Teixeira Pinto – CRF/SP 42654

ANVISA Registration number: 80829430003

SAC 0800-709 6379

If you would like a printed version of this material, at no additional cost, please contact our SAC.

Copyright © 2023 Merz North America, Inc. RADIESSE is a registered trademark of Merz North America, Inc.

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Copyright © 2023 Merz North America, Inc. RADIESSE é uma marca registrada da Merz North America, Inc.