# RADIESSE® Duo

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---

### DESCRIPTION

RADIESSE® Duo injectable implant is a subdermal, semi-solid, cohesive and resorbable implant, non-pyrogenic, latex free, sterile with steam. The principal component is synthetic calcium hydroxyapatite, a biocompatible material. The semi-solid nature of the implant originates from the suspension of calcium hydroxyapatite in a gel-based vehicle, consisting mainly of water (sterile water for injection USP), sodium carboxymethylcellulose (USP) and glycerin (USP). The CaHA-containing gel is formulated to stimulate an increase in the production of collagen and elastin in the treated region. The elastic and collagen fibers are in the middle of the extracellular matrix (ECM) of the dermis and are responsible for important biomechanical properties of the skin. The result is a long-lasting restoration and volumization, with improvements in laxity and skin tightening in the injected region by collagen biostimulation, though not permanent.

*RADIESSE® Duo injectable implant 0.8cc, 1.5cc and 3.0cc have a particle size range of 25–45 microns. The 0.8cc and 1.5cc presentations should be injected with a 27G needle, and the 3.0cc with a 25G needle, with a standard Luer fitting. The use of needles of caliber less than 27G may increase the incidence of occlusion.*

### PRESENTATIONS

RADIESSE® Duo is supplied in a sterilized filled syringe for single use. Each carton contains one syringe of 0.8cc, 1.5cc or 3.0cc, one instruction for use, 3 traceability labels, two sterile needles, for single use only. The dimensions of the needles are specified on the outer carton.

### COMPOSITION

The RADIESSE® Duo Injectable Implant consists of approximately 30% calcium hydroxyapatite per volume and 70% gel as vehicle, by volume.

Composition: calcium hydroxyapatite, sterile water for injection, glycerin and sodium carboxymethylcellulose, sodium.

### INTENDED USE

The RADIESSE® Duo Injectable Implant is indicated for:

- Plastic and reconstructive surgery, including deep filling of the subdermal soft tissues such as correction of wrinkles and moderate to severe furrows (such as nasolabial folds), contouring and improvement in laxity and skin tightening in the injected region by biostimulation of collagen.
- Restoration and/or correction of the signs of facial fat loss (lipoatrophy) in people with human immunodeficiency virus (HIV).
- Correction of volume loss in the dorsum of the hands, rejuvenation of the hands.

It acts by adding the volume of tissue, restoring or increasing the contours of the skin, depending on the desired amount of correction.

The CaHA gel matrix is formulated to stimulate an increase in the production of collagen in the injected region with improvements in laxity and skin tightening (areas such as facial, abdomen, thighs and brachial zone, upper arms, neck and décolletage, etc.). Collagen and elastin are components of the extracellular matrix (ECM) that confer important biomechanical properties to the skin. The result is a long-lasting restoration and volumization, with improvements in laxity and skin tightening in the injected region by collagen biostimulation, though not permanent.

### CONTRAINDICATIONS

RADIESSE® Duo Injectable Implant is contraindicated for patients with:

- severe allergies manifested by a history of anaphylaxis, or history or presence of multiple severe allergies.
- if the patient has active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- in patients with Known hypersensitivity to any components of the formulation.
- in patients prone to developing inflammatory skin conditions or those patients with a tendency for developing hypertrophic scars and keloids.
- implant in the epidermis or use as a skin replacement. Implantation into the epidermis or superficial dermis could lead to complications such as fistula formation, infections, extrusions, nodule formation and induration.

- is not intended to be used for the correction of glabellar folds. A higher incidence of localized necrosis has been associated with glabellar injection. Complications associated with injectables indicate that forceful injection into superficial dermal vessels of the glabellar can result in vascular occlusion.
- in the presence of foreign bodies such as liquid silicone or other particulate materials.
- in areas where there is inadequate coverage of healthy, well vascularized tissue.
- in patients with systemic disorders which cause poor wound healing or will lead to tissue deterioration over the implant.
- patients with bleeding disorders.

## WARNINGS

- *Special attention has been taken in order to avoid the implant into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention.*
- Use of RADIESSE® Duo injectable implant in any person with active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- Do not overcorrect (overfill) a contour deficiency because the depression should gradually improve within several weeks as the treatment effect of RADIESSE® Duo injectable implant occurs. Refer to the individualization of treatment sections for details.
- The safety and effectiveness for use in the lips has not been established. There have been published reports of nodules associated with the use of RADIESSE® Duo injectable implant injected into the lips.
- Special care should be taken to avoid injection into veins or tendons in the hand. Injection into tendons may weaken tendons and cause tendon rupture. Injection into veins may cause embolization or thrombosis.
- Injection procedure reactions have been observed consisting mainly of short-term (i.e., < 7 days) bruising, redness and swelling. Refer to the Adverse Events sections for details.
- Injection into the dorsum of the hand may cause adverse events that last for more than 14 days. Refer to adverse events sections for details.
- Injection in the dorsum of the hand may result in temporary difficulty performing activities (48% of study patients reported this adverse event). Fitzpatrick Skin Types IV-VI may have an increased risk in difficulty performing activities (68% of Fitzpatrick Skin Types IV-VI reported this event).
- RADIESSE® Duo may cause nodules, bumps or lumps in the dorsum of the hand (12% reported this event) and can last up to a 1 year.
- Injection into patients with very severe loss of fatty tissue with marked visibility of veins and tendons has not been studied. The safety and effectiveness in this patient population has not been established.
- Volumes over 3cc of RADIESSE® Duo per hand in a treatment session have not been studied. Increased bruising is associated with higher volume injection. Re-treatment with RADIESSE® Duo of volumes greater than approximately 1.6cc per hand in a treatment session can result in increased adverse events (redness, pain, swelling, and difficulty performing activities).

## PRECAUTIONS

- In order to minimize the risks of potential complications, this product should only be used by healthcare practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. Such practitioners should familiarize themselves with the product, the educational materials of the product and its instructions of use. Healthcare practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Universal precautions must be observed when there is a potential for contact with patient body fluids. The injection session must be conducted with aseptic technique.
- Packaged for single patient use. Do not resterilize. Do not use if package is opened or damaged. Do not use if the syringe end cap or syringe plunger is not in place.
- The calcium hydroxylapatite (CaHA) particles of RADIESSE® Duo injectable implant are radiopaque and are clearly visible on CT Scans and may be visible in standard, plain radiography. In a radiographic study of 58 faces, there was no indication that RADIESSE® Duo injectable implant potentially masked abnormal tissues or being interpreted as tumors in CT Scans. Patients need to be informed of the radiopaque nature of RADIESSE® Duo injectable implant, so that they can inform their primary care health professionals as well as radiologists. Imaging studies have not been performed in the hand. It is presently unknown if RADIESSE® Duo could mask a hand injury on imaging studies.
- Safety of RADIESSE® Duo injectable implant beyond 3 years has not been investigated in clinical trials.

- The safety of RADIESSE® Duo injectable implant in patients with increased susceptibility to keloid formation and hypertrophic scarring has not been studied.
- As with all transcutaneous procedures, RADIESSE® Duo injectable implant injection carries a risk of infection. Standard precautions associated with injectable materials should be followed.
- Safety of RADIESSE® Duo injectable implant for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established.
- The patient should be informed that he or she should minimize strenuous activity and exposure of the treated area to extensive sun or heat exposure for approximately 24 hours after treatment and until any initial swelling and redness has resolved.
- Patients who are using medications that can prolong bleeding, such as aspirin or warfarin, may, as with any injection, experience increased bruising or bleeding at the injection site.
- Safety and effectiveness in the periorbital area has not been established.
- No studies of interactions of RADIESSE® Duo injectable implant with drugs or other substances or implants have been conducted.
- The safety of RADIESSE® Duo injectable implant with concomitant dermal therapies such as epilation, UV irradiation, or laser, mechanical or chemical peeling procedures has not been evaluated in controlled clinical trials.
- If laser treatment, chemical peeling, or any other procedure based on active dermal response is considered after treatment with RADIESSE® Duo injectable implant, there is a possible risk of eliciting an inflammatory reaction at the implant site. This also applies if RADIESSE® Duo injectable implant is administered before the skin has healed completely after such a procedure.
- Injection of RADIESSE® Duo injectable implant into patients with a history of previous herpetic eruption may be associated with reactivation of the herpes.
- Use of RADIESSE® Duo in the dorsum of the hand may result in significant swelling of the dorsum of the hand. Patients should be instructed to remove jewelry (rings) before treatment and until swelling has resolved to avoid compromise of finger circulation.
- Use of RADIESSE® Duo in the dorsum of the hand in patients with diseases, injuries or disabilities of the hand has not been studied. Care should be used in treating patients with autoimmune disease affecting the hand, hand implants, Dupuytren's contracture, history of hand tumor, vascular malformations, Raynaud's disease and patients at risk for tendon rupture.
- To help avoid needle breakage, do not attempt to straighten a bent needle. Discard it and complete the procedure with a replacement needle.
- After use, treatment syringes and needles may be potential biohazards. Handle accordingly and dispose of in accordance with accepted medical practice and applicable local, state and federal requirements.

## **ADVERSE EVENTS**

Adverse events observed in clinical trials with the RADIESSE® Duo Injection Implant include reactions due to injection, needle.

### WRINKLES AND FOLDS

In a randomized, controlled, multicenter study with RADIESSE® Duo injectable implant to the nasolabial folds treatment, very common adverse events (> 10% of patients) reported in the first 14 days were: edema, erythema, ecchymosis, pain and pruritus. Uncommon adverse events (between 0,1% and 1% of patients) were nodule formation.

In a randomized, multicenter study with RADIESSE® Duo injectable implant to for the treatment of facial lipoatrophy associated to HIV, very common adverse events (> 10% of patients) reported in the firsts 14 days were: edema, ecchymosis, erythema, pain, pruritus and protuberance. Common adverse events (between 1% and 10% of patients) were: hardness, discoloration and headache. The events reported in both studies were generally mild, of short duration, and resolved without clinical sequelae.

It was performed a post marketing study to collect long-term safety information of RADIESSE® Duo injectable implant to assess the effect of multiple injections. There were no reports of long term adverse events in this post approval study. The adverse events monitored in the post-approval study included allergic reaction, ecchymosis, edema, embolization, erosion, erythema, extrusion, granuloma, hematoma, infection, necrosis, needle jamming, nodule, and pain.

### FITZPATRICK SKIN TYPE IV-VI STUDIES

Adverse events reported in the short-term Fitzpatrick Skin Type IV-VI post-approval study (100 patients) are erythema, edema, and ecchymosis.

### DORSUM OF THE HAND

A randomized, masked, multicenter study with RADIESSE® Duo injectable implant to dorsum of the hands injections included as very common adverse events (> 10% of patients) lasting more than 14 days were: swelling and pain. Common adverse events (between 1% and 10% of patients) were: nodules, lumps and difficulty performing activities (hands), redness, bruising and hematoma. All events resolved without sequelae.

#### Adverse Events with Duration Greater Than 14 Days:

Events reported by subjects and/or Health professionals to last for longer than 14 days are listed below. The percentages are the number of subjects that experienced an adverse event for greater than 14 days out of 113 subjects that were treated in the study. All events resolved without sequelae.

- 29% swelling; • 25% pain; • 7% nodules/bumps/lump; • 6% difficulty performing activities; • 6% redness; • 3% bruising; • 1% hematoma.

#### POST MARKETING SURVEILLANCE

The following adverse events were received from post surveillance for the injectable RADIESSE® Duo implant. They were not observed in clinical trials. They are voluntary reports so it is not possible to estimate their frequency or to establish a causal relationship with RADIESSE® Duo: pruritus, itching, angioedema, inflammation, necrosis, granuloma, nodules, induration, erythema, skin discoloration, pustule, hair loss, paresthesia, ptosis, pain, headache, swelling, asymmetry, abscess, herpetic infection including herpes simplex and herpes zoster, hematoma, blanching, blistering, dizziness, flu-like symptoms, Guillain-Barre syndrome, tachypnea, ischemic reaction, lymphoid hyperplasia, nausea, skin pallor, scarring, sensitivity to cold, altered skin texture, vascular compromise, ocular ischemia, diplopia, product displacement, facial muscle paralysis, Bell's palsy.

The following interventions have been reported: antibiotics, anti-inflammatories, corticosteroids, anti-histamines, analgesics, massage, warm compress, excision, drainage, and surgery. This information does not constitute and is not intended to be medical advice, a recommendation on how to treat an adverse event or an exhaustive list of possible interventions. Health professionals should evaluate each case on an individual basis, and independently determine, based on their professional experience, what treatment(s) are appropriate, if any, for their patients.

#### INDIVIDUALIZATION OF TREATMENT

Before treatment, the patient's suitability for the treatment and the patient's need for pain relief should be assessed. The outcome of treatment with RADIESSE® Duo injectable implant will vary between patients. In some instances, additional treatments may be necessary depending on the size of the defect and the needs of the patient.

#### DIRECTIONS FOR USE

##### GENERAL

The following is required for the percutaneous injection procedure:

- RADIESSE® Duo injectable implant syringe(s)
  - 25 gauge OD - 27 gauge ID needle(s) with Luer lock fittings
1. Prepare patient for percutaneous injection using standard methods. The treatment injection site should be marked with a surgical marker and prepared with a suitable antiseptic for each patient. Local or topical anesthesia at the injection site or sedation should be used at the discretion of the Health professional.
  2. Prepare the syringes of RADIESSE® Duo injectable implant and the injection needle(s) before the percutaneous injection. A new injection needle may be used for each syringe, or the same injection needle may be connected to each new syringe for each patient.
  3. Remove foil pouch from the carton. Open the foil pouch by tearing at the notches (marked 1 and 2) and remove the syringe from the foil pouch. *There is a small amount of moisture normally present inside the foil pouch for sterilization purposes; this is not an indication of a defective product.*
  4. Peel or twist apart the needle packaging to expose the hub.
  5. Remove the Luer syringe cap from the distal end of the syringe prior to attaching the needle. The syringe of RADIESSE® Duo injectable implant can then be twisted onto the Luer lock fitting of the needle taking care not to contaminate the needle. Discard needle package. The needle must be tightened securely to the syringe and prepared with RADIESSE® Duo injectable implant. If excess implant is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until RADIESSE® Duo injectable implant extrudes from the end of the needle. If leakage is noted at the Luer fitting, it may be necessary to tighten the needle, or to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the needle.
  6. Locate the initial site for the implant. Scar tissue and cartilage may be difficult or impossible to treat. Avoid passing through these tissue types when advancing the injection needle during the procedure.
  7. The amount injected will vary depending on the site and extent of the restoration or augmentation desired. Should be injected subdermally. For rejuvenation of the hands, RADIESSE® Duo injectable implant must be injected at the areolar plane between the subcutaneous layer and superficial fascia in the dorsum of the hand.
  8. Use a 1:1 correction factor. No overcorrection is needed.

9. Insert needle with bevel down at approximately a 30° angle to the skin. Needle should slide under the dermis to the point you wish to begin the injection. This should be easily palpable with the non-dominant hand.
10. If significant resistance is encountered when pushing the plunger, the injection needle may be moved slightly to allow easier placement of the material or it may be necessary to change the injection needle. Needle jamming occurred in the nasolabial folds clinical study. Needle jams are more likely with use of needles smaller than 27 gauge ID.
11. Advance the needle into the subdermis to the starting location. Carefully push the plunger of the RADIESSE® Duo injectable implant syringe to start the injection and slowly inject the implant material while withdrawing the needle in linear threads while withdrawing the needle. Continue placing additional lines of material until the desired level of correction is achieved. For hands indication, inject in small boluses, (0.2 – 0.5cc/bolus) between the central tendons in the dorsum of the hand. Skin tenting should be performed to separate the skin from vascular and tendinous structures and to distribute the product for later massaging.
12. Apply slow continuous even pressure to the syringe plunger to inject the implant as you withdraw the needle. The implant material should be completely surrounded by soft tissue without leaving globular deposits. The injected area may be massaged as needed to achieve even distribution of the implant. For hands indication, apply slow continuous uniform pressure on the plunger of the syringe to inject the implant to deposit the bolus and distribute the product for later massaging.
13. Use once and discard in accordance with local safety standards.


#### Technique for Mixing RADIESSE® Duo injectable implant and 2% Lidocaine HCl

CAUTION: Do not use the RADIESSE® Duo injectable implant and 2% lidocaine mixture later than 2 hours after mixing.

CAUTION: The assembled components are intended for one-time use only.

#### Component Assembly and Mixing Instructions


1. Assemble the components and perform the mixing using sterile technique (see Figure 1).


**Figure 1:** Left to right: Female-to-female luer lock connector, RADIESSE® Duo syringe, mixing syringe, sterile 27 gauge, 0.5" or 25 gauge, 0.5" needle.

**Neither the lidocaine nor the 3.0cc mixing sterile needle nor the Luer lock connector is supplied with this product.**

2. Draw the lidocaine into a sterile polypropylene mixing syringe fitted with a sterile 27 gauge, 0.5" or 25 gauge, 0.5" needle.
3. Tap the mixing syringe, containing lidocaine and depress its push rod to remove all excess air.
4. Remove the sterile 27 gauge, 0.5" or 25 gauge, 0.5" needle.
5. Firmly connect the mixing syringe to the RADIESSE® Duo syringe using the female-to-female luer lock connector (see Figures 2 and 3).


**Figure 2**

**Figure 3**

6. Mix the lidocaine and RADIESSE® Duo injectable implant by alternately depressing the plungers, first on the mixing syringe and then on the RADIESSE® Duo syringe for ten mixing strokes (each mixing stroke is one complete compression of the mixing syringe plunger followed by one complete compression of the RADIESSE® Duo syringe plunger). Plungers are compressed firmly and quickly, at about two compressions per second.

**Figure 4**

7. After mixing, remove the mixing syringe and the female-to-female luer lock connector and discard.
8. Fit the syringe containing the lidocaine and RADIESSE® Duo mixture with an injection needle.
9. Proceed with the injection of the RADIESSE® Duo injectable implant.

The clinical study was conducted by mixing 0.2cc of 2% lidocaine with 1.3cc of RADIESSE® Duo injectable implant in the 3.0cc syringe. The table below provides the ratio of 2% lidocaine to be mixed with the various syringe volumes of RADIESSE® Duo injectable implant. These ratios result in the same concentration of 2% lidocaine (w/v%) in RADIESSE® Duo injectable implant that was mixed in the clinical study after accounting for the dead space in the RADIESSE® Duo and 3.0cc mixing syringes.

Table. LIDOCAINE CONCENTRATION

RADIESSE® Duo (cc)	2% de lidocaine (cc)	Resulting Lidocaine Concentration (w/v%)
0,8	0,11	0,31% - 0,32%
1,5	0,26	0,31% - 0,32%
3,0	0,45	0,32% - 0,34%

**Note:** Health professionals should evaluate each case, individually, and independently determine, based on their professional experience, the appropriate treatment (s), if any, for their patients. Thus, the dilution of CaHA to the treatment may change according to the area treated, the degree of laxity and the thickness of the skin.

#### INSTRUCTIONS FOR TRACEABILITY LABELS

A set of labels is supplied inside the package. The labels on the aluminum packaging are used to ensure product identification and traceability. These labels should be affixed as follows: one in the medical record, another in the document to be delivered to the patient, and the third in the control of the clinic / hospital.

#### STORAGE

RADIESSE® Duo injectable implant should be stored at a controlled room temperature between 15°C and 32°C (59°F and 90°F). Do not use if the expiration date has been exceeded. The expiration date is printed on the product label.

#### PATIENT COUNSELING INFORMATION

For about 24 hours after injection, the following should be avoided in the treated area;

- Significant movements or massages in the treated area
- Applying makeup
- Extensive exposure to sun or heat

In addition, the health professional's instructions regarding care of the injection site should be followed and may include the use of cold packs and / or medications to reduce swelling and bruising.

#### SUPPLIER

The RADIESSE® Duo Injectable Implant is supplied sterile and non-pyrogenic in a syringe packaged in an aluminum foil pack, placed inside a carton for easy storage. Each unit includes a 1.5cc, 0.8cc or 3.0cc, filled syringe of RADIESSE® Duo Injectable Implant. The degree of accuracy of syringe graduations is  $\pm 0.025$ cc. Do not use the product if the package and / or syringe is damaged or if the protective cap or plunger of the syringe is not intact.

**The contents of the syringe are intended for only one patient in a single treatment and cannot be resterilized. Reuse can compromise the functional properties of the device and / or lead to device failure. Reuse can also create a risk of contamination of the device and / or cause infection or cross-infection including, but not limited to, transmission of infectious diseases and transfer of blood between patients. All this can, in turn, lead to injury, illness or death of the patient.**

## DISPOSAL

Used and partially used syringes and injection needles could be biohazardous and should be handled and disposed of in accordance with facility medical practices and local, state or federal regulations.

## WARRANTY

Merz North America, Inc. warrants that reasonable care has been exercised in the design and manufacture of this product.

**THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ITS PARTICULAR PURPOSE.**

Handling and storage of this product, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Merz North America, Inc.'s control directly affect the product and the results obtained from its use. Merz North America, Inc.'s obligation under this warranty is limited to the replacement of this product and Merz North America, Inc. shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly, arising from the use of this product. Merz North America, Inc. neither assumes, nor authorizes any person to assume for Merz North America, Inc., any other or additional liability or responsibility in connection with this product.

## SYMBOLS

	Caution.		Consult instructions for use.
	Do not use if package is damaged.		Single use product. Do not re-use. Forbidden reprocess.
	Not Made with Natural Rubber Latex.		Use by date.
	Date of manufacture.		Manufacturer.
	Lot.		Sterile. Sterile using steam.
	32°C (90°F) 15°C (59°F)	Temperature limit of storage: 15°C – 32°C	

## Manufacturer

Merz North America, Inc.  
4133 Courtney St., Suite 10 - Franksville, WI 53126 USA  
Country of Origin: United States of America.

## Registered/ Imported and Distributed by: Merz Farmacêutica Comercial Ltda.

Headquarters: CNPJ: 11.681.446/0001-45 Avenida Doutor Cardoso de Melo nº 1184 Conj.131 e 132, Vila Olímpia, São Paulo/SP CEP 04548-004  
Branches: CNPJ: 11.681.446/0003-07  
CNPJ: 11.681.446/0005-79

Technical Responsible: Luciana S. M. Franco Silva - CRF/SP 68.932  
ANVISA Registration number: 80829430004  
SAC 0800-709 6379

If you would like a printed version of this material, at no additional cost, please contact our SAC.

Copyright © 2023 Merz North America, Inc. RADIESSE is a registered trademark of Merz North America, Inc.

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Copyright © 2023 Merz North America, Inc. RADIESSE é uma marca registrada da Merz North America, Inc.